


Changing Objectives (u880/i880)

Setup

1. Move nosepiece to an empty position or unscrew one of the current objectives (if all positions are full). Make sure you put the objective taken out in the objective cabinet.
2. Screw in your desired objective. Please handle the objective with care.
3. In Zen, go to the *Maintain* tab.
4. Click on *Objectives* tab.
5. Check the nosepiece position of your objective, then under Detection Objective, click on the objective icon in that particular position.
6. Under *Change Objective*, click on *Favorite Objectives*. You will see a list of objectives.
7. Choose your objective (name, magnification, NA, immersion), then on the menu on the left, click on *Select Objective*.
8. Click *Write to Hardware*.
9. When you go back to the Locate tab, your objective name should be displayed correctly.



Objectives tool for Adding/Changing Objectives

When done

1. Clean objective with lens paper
2. Leave objective in its current position in the nosepiece.