


PALM User Protocols

Chromosome Preparation



PALM Laboratories
A unique service of Carl Zeiss



We make it visible.

PALM User Protocols

Chromosome preparation


1	Introduction	3
2	Consumables and components to work with chromosomes	4
2.1	PALM MembraneSlides	4
2.2	Collectors	5
2.2.1	CapMover	5
2.2.2	RoboMover	5
2.3	Collection devices	6
2.3.1	PALM AdhesiveCap	6
2.3.2	Other microfuge tubes	6
2.3.3	AmpliGrid AG480F	6
3	Chromosome Preparation	7
3.1	Chromosome Preparation from peripheral blood cells	7
3.2	Chromosome Preparation from cell lines (EJ28, HepG2, HeLa)	9
3.3	“Chromosome dropping”	11
4	Chromosome staining: GTG staining (Giemsa Staining)	12
5	Non-contact Laser Capture Microdissection (LCM)	13
5.1	Collection procedures	13
5.1.1	“Dry” collection into PALM AdhesiveCap	13
5.1.2	“Wet” collection into other microfuge tubes	13
5.1.3	Collection onto AmpliGrid AG480F	14
6	Literature	15

Chromosome preparation

1 Introduction

Laser capture microdissection (LCM) is a process to collect pure cell populations for subsequent DNA and RNA extraction. Although laser microdissection is generally used to isolate specific cells from fixed tissue sections, it has been also effective for isolation of individual chromosomes. Preparation of chromosome paints (Kubičková et al. 2002), FISH hybridization (Langer et al. 2005) and DOP-PCR (Hobza et al. 2004) combined with LCM are already described.

One challenge of microdissecting chromosomes is the minimal diameter of samples collected. Non-contact LCM transfers the dissected pure material directly into a collection device. This enables the fast procurement of a homogeneous specimen of just 0.5 μm in diameter without intrusion into the adjacent area.



I.



II.



III.



IV.

Non-contact Laser Capture Microdissection of a metaphase chromosome

- I. human metaphase chromosomes fixed onto PEN membrane (objective: 100x oil)
- II. cutting of selected chromosome (objective: 100x oil)
- III. non-contact LCM of selected chromosome by membrane lifting (objective: 100x oil)
- IV. isolated region after lifting

PALM User Protocols

Chromosome preparation

2 Consumables and components to work with chromosomes

2.1 PALM MembraneSlides

MembraneSlide 0.17 PEN	Order No. 415101-4401-500
MembraneSlide 0.17 PET	Order No. 415101-4401-052
MembraneSlide 1.0 PEN	Order No. 415101-4401-000
MembraneSlide 1.0 PET	Order No. 415101-4401-050
FrameSlide PET	Order no. 415101-4401-055

MembraneSlides are glass slides covered with a membrane on one side. This membrane is easily cut together with the sample and acts as a stabilizing backbone during lifting. Therefore even large areas can be lifted by a single laser pulse without affecting the morphological integrity.

Due to the short working distance of routine 100x (oil) and 63x (oil) magnifying objectives only 0.17 mm thin membrane (PEN, PET) covered glass slides can be used in combination with these objectives!

Note: Thin (0.17 mm) MembraneSlides are not heat-resistant.

An alternative strategy is to use long distance objectives for higher magnifications (20xLD, 40xLD or 63xLD, **not available for 100x!**). There you have to adapt the working distance to the different glass slides by moving the correction collar on the objective to the right position (see figure 2).




Figure 2

left image:	MembraneSlide	1.0 mm	->	recommended setting	= 1
right image:	MembraneSlide	0.17 mm	->	recommended setting	= dot
without image:	FrameSlide		->	recommended setting	= between dot and 0

Chromosome preparation

Carl Zeiss MicroImaging (CZMI) offers 1 mm and 0.17 mm slides covered with polyethylene naphthalate (PEN)-membrane (MembraneSlide PEN). This PEN-membrane (2 µm thickness) is highly absorptive in the UV-A range, which facilitates laser cutting. The PEN-membrane can be used for all kind of applications.

In addition to MembraneSlide PEN, CZMI offers polyethylene teraphthalate-membrane covered (MembraneSlide PET) slides for special processes. The PET-membrane is 1.4 µm thick and highly absorptive in the UV-A range, too. The PET-membrane does not exhibit any autofluorescence itself and is therefore dedicated excellently for fluorescence applications.

It is also suitable to perform laser ablation of unwanted specimen. Any downstream applications are feasible.

Alternatively the PET-membrane can be attached to a metal frame (FrameSlide PET). In fluorescence applications (FISH) even weak signals can be detected due to low signal to noise ratio. The frame structure of FrameSlide PET enables immunohistological staining with minimal volume (< 1 ml) and is resistant to microwave treatment. The special bonding is inert and adapted to heat treatment (up to 95°C) so that the membrane does not ruffle during the heating process.

2.2 Collectors

2.2.1 CapMover

PALM CapMover is designed to be equipped with a SingleTube Collector 200/500 or SingleCap Collector 200/500 or TubeCollector 2x200/2x500 for 200 µl and 500 µl tubes/caps.

2.2.2 RoboMover

PALM RoboMover is a device for automated harvesting and sorting of different kinds of microdissected specimen in a higher throughput process. The non-contact LCM process can be operated automatically or manually. The RoboMover can be reconfigured with different collectors, e.g. SingleTube Collector, TubeCollector 12x8, CapturePlate, Collector 96 or SlideCollector 48.

PALM User Protocols

Chromosome preparation

2.3 Collection devices

2.3.1 AdhesiveCap

AdhesiveCap 500 clear
AdhesiveCap 200 clear

Order No. 415101-4400-255
Order No. 415101-4400-245

The intention of AdhesiveCaps is to allow non-contact LCM without applying any capturing liquid into the caps.

Beside the quick relocation of the lifted samples in the cap due to instant immobilization there is no risk of evaporation and crystal forming of buffers during extended specimen harvesting. For more details and handling, please see also AdhesiveCap product information.


Note: Carl Zeiss recommends AdhesiveCaps as a collection device for most experiments.

2.3.2 Other microfuge tubes

In case you use commercially available plasticware (e.g. ABgene #AB-0350; 0.5 ml tubes) without adhesive filling, a capturing liquid is necessary to attach the microdissected samples.

2.3.3 AmpliGrid AG480F

Using the SlideCollector 48 in conjunction with AmpliGrid technology from Advalytix enables working in a higher throughput LCM. The AmpliGrid technology allows DNA amplification and sequencing in an extremely low volume (1 μ l) directly on chip.



(Please see brochure: On-chip Single Cell Analysis for PALM MicroBeam; to order at ZEISS with brochure no. 60-3-0002/e)

PALM User Protocols

Chromosome preparation

3 Chromosome Preparation

3.1 Chromosome Preparation from peripheral blood cells

Reagents:

RPMI Medium	(Biochrom # FG1385)
Fetal Bovine Serum	(Biochrom # S0113)
Penicillin/Streptomycin: 10000µg/ml	(Biochrom # A2213)
Phytohemagglutinin	(Biochrom # M5030)
Colcemid	(Biochrom # L6231)
Hypotonic Solution: 0.075 M KCl	
Fixation Mix: Methanol/Glacial Acetic Acid (3:1)	

Labware:

Falcon tubes (15 ml)
Glass pipettes (230 mm) and sucker
MembraneSlides (0.17 mm or 1 mm, PET or PEN; FrameSlide PET)
Hettich Centrifuge Universal 32R

Prearrangements:

preheat 0.075 M KCl at 37°C
freshly prepare Methanol/Glacial Acetic Acid (3:1) and store at -20°C
put on ice shortly before use

Protocol:

- set up blood cultures (sterile flask):

RPMI Medium	8 ml
Fetal Bovine Serum	2 ml
Phytohemagglutinin	0.2 ml
Penicillin	0.05 ml
heparinized blood	8 drops ~ 0.3 ml (adult) 5-6 drops ~ 0.2 ml (newborn)

- incubate at 37°C for 72 hours

PALM User Protocols

Chromosome preparation

Colcemid treatment:

sway the flask add Colcemid	100 μ l
--------------------------------	-------------

- incubate at 37°C for 30 min

Note: Prolonged Colcemid treatment may result in shrunken chromosomes.

- after colcemid incubation transfer the suspension into a 10-15 ml tube
- centrifuge at 1000 rpm for 10 min

Metaphase preparation:

remove supernatant resuspend pellet in 6 ml KCl (0.075M), by adding drop by drop during careful mixing at room temperature
--

- incubate at 37°C for 10-20 min
- spin down at 1000 rpm for 10 min
- discard supernatant (**Caution:** pellet is somewhat wobbly)
- resuspend pellet cautiously in some drops of ice-cold Fixation Mix (see Reagents) (stored on ice during the next steps)
- slowly fill up with Fixation Mix to 6 ml

Note: It is very important to resuspend the pellet carefully. This means to pipette the suspension carefully up and down or to mix it by gently knocking at the tube.

- repeat the centrifugation step and resuspend till supernatant gets clear after the last centrifugation step
- resuspend the pellet in about 2-5 ml Fixation Mix

Note: Adjust volume for pellet resuspension to the wanted density of cells and metaphase spreads on the slides. If a high amount of cells is present in the suspension use more Fixation Mix to thin it out. More slides can be prepared this way.

PALM User Protocols

Chromosome preparation

3.2 Chromosome Preparation from cell lines (EJ28, HepG2, HeLa)

Depending on the cell type, adequate media, antibiotics etc. have to be used for the cultivation.

Please, see PALM protocols:

*“Cell Culture, non-contact Laser Capture Microdissection and recultivation” or
“Live Cells and Molecular Analysis”*

Reagents:

Hanks' medium (with Ca ⁺⁺ , Mg ⁺⁺)	(Biochrom # L 2035)
Trypsin / EDTA	(Biochrom # L 2143)
Colcemid	(Biochrom # L 6231)
Hypotonic Solution: Hanks'/Aqua dest. (1:3)	
Fixation Mix: Methanol/Glacial Acetic Acid (3:1)	

Labware:

Falcon tubes (15 ml)
Glass pipettes (230 mm) and sucker
MembraneSlides (0.17 mm or 1 mm; PET or PEN)
Hettich Centrifuge Universal 32R

Prearrangements:

preheat cell culture medium, Hanks' and Trypsin/EDTA at least at room temperature
preheat hypotonic solution at 37°C
freshly prepare Methanol/Glacial Acetic Acid (3:1) and store at -20°C;
put on ice shortly before use

Protocol:

- set up cell cultures (sterile flask):

RPMI Medium	9 ml
Fetal Bovine Serum	0.9 ml
Penicillin	0.09 ml
cell suspension	20 µl (depending on the cell type)

- incubate at 37°C for 72 hours (depending on the cell type)

PALM User Protocols

Chromosome preparation

Colcemid treatment:

4-16 hours before culture stop add Colcemid	100-200 μ l
--	-----------------

Note: Prolonged Colcemid treatment may result in shrunken chromosomes

Trypsinization:

remove medium	
rinse with Hanks' solution	5 ml
remove the solution	
add Trypsin/EDTA	3 ml

- incubate at 37°C for 5 minutes

Note: If cells are hard to trypsinize, knock gently at the flask

to stop trypsinization: add medium transfer the suspension into a Falcon tube	8 ml
---	------

- spin down at 1000 rpm for 10 minutes

Metaphase preparation:

remove supernatant resuspend pellet in 8 ml KCl (0.075M), containing Hanks'/Aqua dest. (1:3) by adding it drop by drop
--

- incubate at 37°C for 30-40 min
- spin down at 800-1000 rpm for 10 min
- discard supernatant (**Caution:** pellet is somewhat wobbly)
- resuspend pellet carefully in some drops of ice-cold Fixation Mix (stored on ice during the next steps)
- slowly fill up with Fixation Mix to 4 ml

Note: It is very important to resuspend the pellet carefully. This means to pipette the suspension carefully up and down or to mix it by gently knocking at the tube.

PALM User Protocols

Chromosome preparation

- repeat the centrifugation step and resuspend till supernatant gets clear after the last centrifugation step
- resuspend the pellet in about 2-5 ml of Fixation Mix

Note: Adjust volume for pellet resuspension on the wanted density of cells and metaphase spreads on the slides. If a high amount of cells is present in the suspension use more Fixation Mix to thin it out. More slides can be prepared this way.

3.3 “Chromosome Dropping”

- dip a plain slide into 100% ethanol, rinse in distilled water and let air dry
- store the so prepared slides at -20°C until starting to drop the chromosome suspension
- take out a slide immediately before dropping the chromosomes
(Defrosting causes a wet slide and makes spreading much easier.)

Note: If using a routine 1 mm glass slide underneath the 0.17 mm MembraneSlide for “stabilization”, the chromosomes will not bounce during dropping, thus improving chromosome spreading.

- put about 2-4 drops of fixed cell suspension from a distance of about 20 cm onto an ice-cold and wet MembraneSlide
- cover a heating block with wet tissues
- let the slides dry on the heating block for 1 hour at 37°C
This so called ‘chromosome aging’ has direct influence on the quality of the metaphases.

PALM User Protocols

Chromosome preparation

4 Chromosome staining: GTG staining (Giemsa staining)

Reagents:

Giemsa Stain Solution	(Sigma Aldrich # GS 500)
Di-sodium hydrogen phosphate solution (1/15 mol/l)	(Merck # 1.065871)
Potassium dihydrogen phosphate solution (1/15 mol/l)	(Merck # 1.04875)
Trypsin (1:250)	(PAA # L11-002)
70% Ethanol	

Prearrangements (prepare all solutions freshly):

Giemsa Stain:

Di-sodium hydrogen phosphate solution	50 ml
Potassium dihydrogen phosphate solution	50 ml
Giemsa Stain Solution	8 ml

Trypsin solution:

PBS (1x)	100 ml
Trypsin (1x)	3 ml

Note: Use only freshly prepared trypsin!

Protocol:

- soak slide into Trypsin 5 seconds
- dip slide into 96 % Ethanol
- dip slide into Aqua dest
- soak slide into Giemsa stain for 4 minutes
- dip slide into Aqua dest
- let the slide air dry

Note: The quality of chromosome spreads is highly depending on the time period of trypsinization and the time period of previous chromosome aging.

Chromosome preparation


5 Non-contact Laser Capture Microdissection (LCM)

5.1 Collection procedures

Please, additionally have a look into the PALM MicroBeam user manual.

5.1.1 "Dry" collection into AdhesiveCap

- after LCM into AdhesiveCap add a buffer of your own choice and incubate "upside down" for 2-18 hours depending on chromosome preparation



- subsequently centrifuge the liquid into the tip of the tube and then apply the routine procedure

5.1.2 "Wet" collection into other microfuge tubes

When using "unfilled" routine microfuge tubes it is necessary to fill a liquid into the cap prior to LCM. This improves the adhesion of the captured chromosomes. Igepal CA-630 in the Capturing Buffer (see next page) allows to smear out a small amount of liquid in the whole cap area. Alternatively the chromosomes can be lifted directly and precisely in a 3 μ l droplet of buffer of your own choice.

Note: All kinds of aqueous solutions will run dry after a while.

Chromosome preparation

Prearrangements:

Capturing Buffer:

0.5 M EDTA pH 8.0	20 μ l
1 M Tris pH 8.0	200 μ l
Igepal CA-630 (SIGMA #I-3021)	50 μ l
Proteinase K 20 mg/ml (Qiagen #19131)	100 μ l
ddH ₂ O	fill up to 10 ml

Note: Always prepare a fresh mixture of Capturing Buffer and Proteinase K.

Protocol:


- take an autoclaved cap or tube
- pipette 3-15 μ l of Capturing Buffer or DNase-free water in the middle of the cap
- check the right position of the correction collar (see page 4)
- put the cap/tube into the collector and position it above the selected chromosome
- perform non-contact LCM of selected chromosome or chromosome parts
- put the cap onto the tube
- centrifuge the tube at full speed for 5 minutes (Hettich Centrifuge, Universal 32R)
- add 10-15 μ l Capturing Buffer containing Proteinase K onto the chromosomes (which are now located in the tip of the tube)
- vortex gently, digest for 2 - 18 hours at 55 °C followed by a final heating step at 90 °C for 10 min to inactivate Proteinase K
- if not going on immediately, store the samples in a refrigerator at 4 °C

Note: The time necessary for complete digestion depends on the type and the number of isolated chromosomes.

- after the Proteinase K digest and the inactivation step you may start with your routine downstream application.

5.1.3 “Wet” collection onto AmpliGrid AG480F

A preloading of the 48 ReactionSites of the AmpliGrid with 1 μ l liquid (e.g., 1% Glycerol in water) enables prolongation of the working time and increases adhesion of the captured samples. The LCM process of 48 ReactionSites can be operated automatically and is controlled by PALM RoboSoftware.



Chromosome preparation

6 Literature

PALM Systems

Aubele MM, Cummings MC, Mattis AE, Zitzelsberger HF, Walch AK, Kremer M, Höfler H and Werner M

Accumulation of chromosomal imbalances from intraductal proliferative lesions to adjacent in situ and invasive ductal breast cancer.

Diagn Mol Pathol, 2000, 9,1: 14-19

Aubele M, Cummings M, Walch A, Zitzelsberger H, Nahrig J, Höfler H and Werner M

Heterogeneous chromosomal aberrations in intraductal breast lesions adjacent to invasive carcinoma.

Anal Cell Pathol, 2000, 20,1: 17-24

Aubele M, Mattis A, Zitzelsberger H, Walch A, Kremer M, Welzl G, Höfler H and Werner M

Extensive ductal carcinoma in situ with small foci of invasive ductal carcinoma: evidence of genetic resemblance by CGH.

Int J Cancer, 2000, 85: 82-86

Fominaya A, Linares C, Loarce Y and Ferrer E

Microdissection and microcloning of plant chromosomes.

Cytogenet Genome Res, 2005, 109: 8-14

Fuková I, Traut W, Vitková M, Nguyen P, Kubičková S and Marec F

Probing the W chromosome of the codling moth, *Cydia pomonella*, with sequences from microdissected sex chromatin.

Chromosoma, 2006, 116,2: 135-145

Gotter AL, Nimmakayalu MA, Jalali GR, Hacker AM, Vorstman J, Conforto Duffy D, Medne L and Emanuel BS

A palindrome-driven complex rearrangement of 22q11.2 and 8q24.1 elucidated using novel technologies.

Genome Res, 2007, 174: 470-481

Hartmann A, Schlake G, Zaak D, Hungerhuber E, Hofstetter A, Hofstaedter F and Knuechel R

Occurrence of chromosome 9 and p53 alterations in multifocal dysplasia and carcinoma In Situ of human urinary bladder.

Cancer Res, 2002, 62: 809-818

Chromosome preparation

Hartmann A, Rösner U, Schlake G, Dietmaier W, Zaak D, Hofstaedter F and Knuechel R

Clonality and genetic divergence in multifocal low-grade superficial urothelial carcinoma as determined by chromosome 9 and p53 deletion analysis.

Lab Invest, 2000, 80,5: 709-718

Koene GPJA, Arts-Hilkes YHA, van der Ven KJW, Rozemuller EH, Slootweg PJ, de Weger RA and Tilanus MGJ

High level of chromosome 15 aneuploidy in head and neck squamous cell carcinoma lesions identified by FISH analysis: limited value of β 2-microglobulin LOH analysis.

Tissue Antigens, 2004, 64: 452-461

Kubičková S, Cernohorska H, Musilova P and Rubes J

The use of laser microdissection for the preparation of chromosome-specific painting probes in farm animals.

Chromosome Res, 2002, 10: 571-577

Langer S, Geigl JB, Gangnus R and Speicher MR

Sequential application of interphase-FISH and CGH to single cells.

Lab Invest, 2005, 85,4: 582-592

Luebke AM, Schlomm T, Gunawan B, Bonkhoff H, Füzesi L and Erbersdobler A

Simultaneous tumour-like, atypical basal cell hyperplasia and acinar adenocarcinoma of the prostate: a comparative morphological and genetic approach.

Virchows Archiv, 2005, 446,3: 338-341

Matsunaga S, Kawano S, Michimoto T, Higashiyama T, Nakao S, Sakai A and Kuroiwa T

Semi-automatic laser beam microdissection of the Y chromosome and analysis of Y chromosome DNA in a dioecious plant, *Silene latifolia*.

Plant Cell Physiol, 1999, 40,1: 60-68

McNally LR, Henk WG, Cooper RK

Laser pressure catapulting followed by B actin gene identification in Japanese quail macrochromosomes and microchromosomes using Teflon-coated coverslip slides.

J Microsc, 2005, 218: 219-224

Obermann EC, Junker K, Stoehr R, Dietmaier W, Zaak D, Schubert J, Hofstaedter F, Knuechel R and Hartmann A

Frequent genetic alterations in flat urothelial hyperplasias and concomitant papillary bladder cancer as detected by CGH, LOH, and FISH analyses.

J Pathol, 2003, 199: 50-57

Chromosome preparation

Pyakurel P, Montag U, Castañós-Vélez E, Kaaya E, Christensson B, Tönnies H, Biberfeld P and Heiden T

CGH of microdissected Kaposi's sarcoma lesions reveals recurrent loss of chromosome Y in early and additional chromosomal changes in late tumour stages.

Aids, 2006, 20: 1805-1812

Schermelleh L, Thalhammer S, Heckl W, Pösl H, Cremer T, Schütze K and Cremer M

Laser microdissection and laser pressure catapulting for the generation of chromosome-specific paint probes.

BioTechniques, 1999, 27: 362-367

Stark RW, Rubio-Sierra FJ, Thalhammer S and Heckl WM

Combined nanomanipulation by atomic force microscopy and UV-laser ablation for chromosomal dissection.

Eur Biophys J, 2003, 32: 33-39

Thalhammer S, Langer S, Speicher MR, Heckl WM and Geigl JB

Generation of chromosome painting probes from single chromosomes by laser microdissection and linker-adaptor PCR.

Chromosome Res, 2004, 12: 337-343

Thalhammer S, Stark RW and Schütze K

Laser microdissection of metaphase chromosomes and characterization by atomic force microscopy.

J Biomed Optics, 1997, 2,1: 115-119

Werner M, Mattis A, Aubele M, Cummings M, Zitzelsberger H, Hutzler P and Höfler H

20q13.2 amplification in intraductal hyperplasia adjacent to in situ and invasive ductal carcinoma of the breast.

Virchows Arch, 1999, 435,5: 469-472

PALM User Protocols

Chromosome preparation

DOP-PCR

Di Berardino D, Vozdova M, Kubičková S, Cernohorska H, Coppola G, Coppola G, Enne G and Rubes J

Sexing river buffalo (*Bubalus bubalis* L.), sheep (*Ovis aries* L.), goat (*Capra hircus* L.), and cattle spermatozoa by double color FISH using bovine (*Bos Taurus* L.) X- and Y-painting probes.

Mol Reprod Dev, 2004, 67: 108-115

Hobza R, Lengerova M, Cernohorska H, Rubes J and Vyskot B

FAST-FISH with laser beam microdissected DOP-PCR probe distinguishes the sex chromosomes of *Silene latifolia*.

Chromosome Res, 2004, 12,3: 245-250

Rubes J, Kubickova S, Musilova P, Amaral EM, Brunner RM and Goldammer T

Assignment of chromosome rearrangements between X chromosomes of human and cattle by laser microdissection and zoo-FISH.

Chromosome Res, 2005, 13,6: 569-574

Telenius H, Pelmeur AH, Tunnacliffe A, Carter NP, Behmel A, Ferguson-Smith MA, Nordenskjold M, Pfragner R and Ponder BA

Cytogenetic analysis by chromosome painting using DOP-PCR amplified flow sorted chromosomes.

Genes Chromosomes Cancer, 1992, 4: 257-263

Telenius H, Carter NP, Bebb CE, Nordenskjold M, Ponder BA, Tunnacliffe A.

Degenerate Oligonucleotide-Primed PCR: General Amplification of Target DNA by a single Degenerate Primer.

Genomics, 1992, 13: 718-72